


# YEAST PROPAGATION



**AGAR SLANT OR  
COTTON BALL**



**100 - 200 ML  
STERILE WORT  
AT AMBIENT TEMP**



**LIQUID YEAST IN  
ALUMINIUM  
BOTTLE**

**1 L THEN GO UP  
TO 5 - 10 L  
STERILE WORT AT  
AMBIENT TEMP**



**RESUSPENDED  
PRESSED YEAST**



**CARLSBERG  
FLASK AT  
AMBIENT TEMP**